

Results Cannot Be Approved Article ID: KB-001

Summary

This article explains how to troubleshoot and resolve situations where test results cannot be approved in AcmeLIMS.

Environment

- AcmeLIMS
- Results Management
- Sample Approval

Symptoms

When attempting to approve a sample, the **Approve** button is unavailable or an error message is displayed.

Cause

Results can only be approved when all prerequisite laboratory processes have been completed. Common causes include:

- Mandatory tests have not been completed.
- Results are still in draft status.
- Required validation checks have not been performed.
- The user does not have approval permissions.

Resolution

1. Open the sample record.
2. Review the sample status.
3. Confirm that all required tests have been completed.
4. Verify that the results have been validated.
5. Check that no mandatory fields are blank.
6. Confirm that your user account has the appropriate approval role.
7. Attempt to approve the results again.

Additional Information

Result approval forms part of the laboratory audit trail and should only be performed by authorised personnel.

Related Articles

- KB-002: Validating Test Results
- KB-005: Managing User Roles and Permissions
- KB-008: Understanding Sample Statuses