

AcmeLIMS Laboratory Information Management System

Operational Qualification Test Script

Aliquot Sample

Operational Qualification Test Script Approval

Your signature below indicates your approval of the scope and content of this test script.

Approved by:

Company	Name/Title	Signature	Date
AcmeLIMS			
(Company Name)			
(Company Name QA)			

Document History

Version	Revision Date	Revision Description	Author
1	05 Jul 2026	New Document	Clare Millerchip, Technical Author, AcmeLIMS

Overview of Tests

Test Stage	Other Operations
Test Environment	The AcmeLIMS Regulated Manufacturing system, Sample Aliquot function, Version _____
Objective of Test	The objective of this test is to verify the functionalities in the Sample Aliquot area of AcmeLIMS .
Description of Sub Tests	1. Open the Aliquot Samples screen 2. Create an Aliquot of a Sample 3. Edit an Aliquoted Sample 4. View the Tests Associated with an Aliquoted Sample
Assumptions and Limitations	The hardware and operating system will not be validated.
Prerequisites	The Installation of AcmeLIMS has been successfully completed.

Test Procedures

1. Open the Aliquot Samples Screen

Test No.	Procedure	Expected Result	Actual Result	Pass/Fail	Initials and Date
1.	In the main Laboratory Process Flow screen: 1. Click the Maintenance button 2. Click the Aliquot button.	The Aliquot Samples screen is displayed and contains the following buttons: • Sample Edit • Tests • Aliquot • Exit	☐ As expected ☐ Other than expected	Pass ☐ Fail ☐	
1.a		There is a box titled Search that contains six graphical controls: • SampleCode (text box) • From (date picker) • To (date picker) • Substance (drop-down list) • Source (drop-down list) • Lot Number (text box – defaults to display an asterisk*) There is a text display under these controls that communicates how many samples are found upon searching in the format Samples found x There are two buttons under these controls: • Reset • Search That can be used with the above graphical controls to search for Aliquot samples.	☐ As expected ☐ Other than expected	Pass ☐ Fail ☐	

AcmeLIMS OQ Test Script – Aliquot Sample

Test No.	Procedure	Expected Result	Actual Result	Pass/Fail	Initials and Date
1.b		There is a box in the top of the screen which has the following columns: • Sample • Batch Code • Substance • Date Received • SubmitterCode • Description These columns can sort the data they contain by ascending or descending order when clicked.	☐ As expected ☐ Other than expected	Pass Fail	☐ ☐
1.c		There is a box in the bottom of the screen that is labelled Associated Aliquots which has the following columns: • SampleCode • Sample Size • Date Created • Time Created • Created By These columns can sort the data they contain by ascending or descending order when clicked.	☐ As expected ☐ Other than expected	Pass Fail	☐ ☐

AcmeLIMS OQ Test Script – Aliquot Sample

Test No.	Comments

2. Create an Aliquot of a Sample

Test No.	Procedure	Expected Result	Actual Result	Pass/Fail	Initials and Date
2.	In the Aliquot Samples screen click the Aliquot button.	The Create Aliquot screen is displayed that contains the following fields: • Parent Sample • Count/Fill Remaining • Units • Number of Aliquots to Produce • Size of Aliquot • Amount Remaining After Aliquoting	☐ As expected ☐ Other than expected	Pass ☐ Fail ☐	
2.a		The Create Aliquot screen also contains two buttons: • OK • Cancel	☐ As expected ☐ Other than expected	Pass ☐ Fail ☐	
2.b	The cursor focus should be in the Number of Aliquots To Produce field. 1. Enter a number in this field for example 3. 2. Click the OK button 3. Click the Sample Label button in the Print Aliquot Labels screen. 4. Close the report window. 5. Click the Sample Label button again. 6. Close the report window. 7. Click the Sample Label button again. The Sample label screen only displays when the OK button is clicked.	• The Print Aliquot Labels screen is displayed after clicking OK with the Windows radio button selected as default. • A report is displayed in a new window. • The report window is closed. • The report for the next aliquot, for example the second, displayed in a new window • The report window is closed. • The report for the next aliquot, for example the third, is displayed in a new window. • The three newly created Aliquoted samples are listed in the Associated Aliquots box at the bottom of the screen.	☐ As expected ☐ Other than expected	Pass ☐ Fail ☐	

AcmeLIMS OQ Test Script – Aliquot Sample

Test No.	Procedure	Expected Result	Actual Result	Pass/Fail		Initials and Date
	8. Close the report window.					

Test No.	Comments

3. Edit an Aliquoted Sample

Test No.	Procedure	Expected Result	Actual Result	Pass/Fail	Initials and Date
3.	In the Aliquot Samples screen click the Sample Edit button.	The Edit Sample screen is displayed. This screen contains the following fields: • BatchCode (text box) • Submitter (selector) • Lot Number (text box) • Date of Sampling/Creation Date (text box) • Sample Container Size (number field and drop down- list) • Sample Description (text box) • Courier (drop – down list) • Due Date (date picker) • SubstanceVersion (text box) • Modification Reason (text box) The BatchCode, Date of Sampling/Creation Date, and Sample Container Size are mandatory fields, as indicated in yellow.	☐ As expected ☐ Other than expected	Pass ☐ Fail ☐	
3.a		The Edit Sample screen also contains the following buttons: • Save • Notes (General) • Cancel	☐ As expected ☐ Other than expected	Pass ☐ Fail ☐	
3.b	1. Change the details in one of the fields in the Edit Sample screen, for example the Submitter name. 2. Click the Save button.	The details for the Aliquot are saved with no errors and displayed correctly in the Aliquot Samples screen.	☐ As expected ☐ Other than expected	Pass ☐ Fail ☐	

AcmeLIMS OQ Test Script – Aliquot Sample

Test No.	Comments

4. View the Tests Associated with an Aliquoted Sample

Test No.	Procedure	Expected Result	Actual Result	Pass/Fail	Initials and Date
4.	In the Aliquot Samples screen click the Tests button.	The Sample Test Management screen is displayed. This screen contains the following buttons: • View • Limits • Associated Files • Add Test • Delete Test • Add Multiple Tests • Allocate Instrument • Allocate Instrument (All Tests) • Exit	☐ As expected ☐ Other than expected	Pass Fail	☐ ☐
4.a		There is a display box labelled Selected Sample and a box labelled Priority above the main box in the screen that displays tests. This lists the current tests that are associated with the aliquoted sample.	☐ As expected ☐ Other than expected	Pass Fail	☐ ☐
4.b		There is a box that contains columns that are labelled the following: • Order • Code • Description • Test Status • Test Position • Test Version • Instrument • Validation Required	☐ As expected ☐ Other than expected	Pass Fail	☐ ☐

AcmeLIMS OQ Test Script – Aliquot Sample

Test No.	Procedure	Expected Result	Actual Result	Pass/Fail		Initials and Date
		Instrument Status These columns can sort the data they contain by ascending or descending order when clicked.				

Test No.	Comments

Test Summary

Test Results Overview		
Test executed by		
Test execution dates	Start Date:	End Date:
Test script execution number:		

Details of Discrepancies

Discrepancy #	Test Step #	Description	Impact	Accepted Y/N
1				
2				
3				
4				
5				
6				
7				
8				
9				
10				
11				
12				
13				
14				

Post Test Script Execution Approval

Your signature below indicates you have reviewed and approve the OQ Test Script execution.

Post OQ Test Script Execution Approval			
Overall Result – cross out as applicable	Pass	Pass with Discrepancy	Fail
Retest required Yes/No			
Additional Comments			

Test Sign Off

Company	Name/Title	Signature	Date
AcmeLIMS			
(Company Name) system owner			
(Company Name) QA			

Appendix – Discrepancy Evidence

Screenshot 1 – Test Step No # [Brief Description]

[Insert Screenshot]

Screenshot 2 - Test Step No # [Brief Description]

[Insert Screenshot]

Screenshot 3 - Test Step No # [Brief Description]

[Insert Screenshot]