

AcmeLIMS Release Notes

Version: 2.5

Release Date: 04 July 2026

Overview

AcmeLIMS Version 2.5 introduces enhancements to result approval workflows, notification management, and sample processing. This release also includes usability improvements and defect fixes designed to improve efficiency and system performance.

New Features

Enhanced Result Approval Workflow

Users can now view validation status directly from the sample record before approving results.

Benefit: Reduces the need to navigate between multiple screens when reviewing samples for approval.

Email Notification Templates

System administrators can configure custom notification templates for:

- Sample status changes
- Approval requests
- Completed testing activities

Benefit: Provides more flexible communication and reduces manual administration.

Improvements

Sample Search Performance

Improved the performance of sample searches when working with large datasets.

Benefit: Reduces response times and improves user productivity.

Audit Trail Visibility

Additional audit information is now displayed when reviewing sample activity.

Benefit: Makes it easier to trace changes and support compliance requirements.

Defect Fixes

Reference	Description
ACM-101	Fixed an issue where email notifications were not sent when sample approvals were completed.
ACM-127	Fixed an issue that prevented users from editing notification rules under certain conditions.
ACM-142	Fixed an issue where completed samples could appear in active worklists.

Known Issues

Users may experience slower performance when generating very large reports containing more than 50,000 records.

A future release will include additional report optimisation.

Upgrade Information

No database changes are required for this release.

Users should ensure all active sessions are closed before the upgrade is performed.

Support

For assistance with this release, contact your system administrator or support team.